

## Ordering Information

### AB-1 (100% Dimethylpolysiloxane)

#### Non Polar

Equivalent

to: BP-1, CP-Sil 5CB, DB-1, HP-1, Ultra-1, Rtx-1,  
SPB-1, ZB-1

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
0153-1008	AB-1, 10m x 0.53mm x 0.88um	Enquire
0153-1026	AB-1, 10m x 0.53mm x 2.65um	Enquire
0153-1050	AB-1, 10m x 0.53mm x 5.0um	Enquire
0120-1203	AB-1, 12m x 0.2mm x 0.33um	Enquire
0125-1501	AB-1, 15m x 0.25mm x 0.10um	Enquire
0125-1502	AB-1, 15m x 0.25mm x 0.25um	Enquire
0125-1510	AB-1, 15m x 0.25mm x 1.0um	Enquire
0132-1501	AB-1, 15m x 0.32mm x 0.11um	Enquire
0132-1502	AB-1, 15m x 0.32mm x 0.25um	Enquire
0132-1510	AB-1, 15m x 0.32mm x 1.00um	Enquire
0153-1501	AB-1, 15m x 0.53mm x 0.15um	Enquire
0153-1502	AB-1, 15m x 0.53mm x 0.25um	Enquire
0153-1505	AB-1, 15m x 0.53mm x 0.50um	Enquire
0153-1510	AB-1, 15m x 0.53mm x 1.00um	Enquire
0153-1515	AB-1, 15m x 0.53mm x 1.50um	Enquire
0153-1530	AB-1, 15m x 0.53mm x 3.00um	Enquire
0153-1550	AB-1, 15m x 0.53mm x 5.00um	Enquire
0120-2501	AB-1, 25m x 0.2mm x 0.11um	Enquire
0120-2503	AB-1, 25m x 0.2mm x 0.33um	Enquire
0120-2505	AB-1, 25m x 0.2mm x 0.50um	Enquire
0132-2501	AB-1, 25m x 0.32mm x 0.17um	Enquire
0132-2505	AB-1, 25m x 0.32mm x 0.52um	Enquire
0132-2510	AB-1, 25m x 0.32mm x 1.05um	Enquire
0125-3001	AB-1, 30m x 0.25mm x 0.1um	Enquire
0125-3002	AB-1, 30m x 0.25mm x 0.25um	Enquire
0125-3005	AB-1, 30m x 0.25mm x 0.50um	Enquire
0125-3010	AB-1, 30m x 0.25mm x 1.00um	Enquire
0132-3001	AB-1, 30m x 0.32mm x 0.10um	Enquire
0132-3002	AB-1, 30m x 0.32mm x 0.25um	Enquire
0132-3005	AB-1, 30m x 0.32mm x 0.5um	Enquire
0132-3010	AB-1, 30m x 0.32mm x 1.0um	Enquire
0132-3030	AB-1, 30m x 0.32mm x 3.0um	Enquire
0132-3040	AB-1, 30m x 0.32mm x 4.00um	Enquire
0132-3050	AB-1, 30m x 0.32mm x 5.00um	Enquire
0153-3001	AB-1, 30m x 0.53mm x 0.11um	Enquire
0153-3002	AB-1, 30m x 0.53mm x 0.25um	Enquire
0153-3005	AB-1, 30m x 0.53mm x 0.50um	Enquire
0153-3008	AB-1, 30m x 0.53mm x 0.88um	Enquire
0153-3010	AB-1, 30m x 0.53mm x 1.00um	Enquire
0153-3015	AB-1, 30m x 0.53mm x 1.50um	Enquire
0153-3026	AB-1, 30m x 0.53mm x 2.65um	Enquire
0153-3030	AB-1, 30m x 0.53mm x 3.00um	Enquire
0153-3050	AB-1, 30m x 0.53mm x 5.00um	Enquire
0120-5001	AB-1, 50m x 0.2mm x 0.11um	Enquire
0120-5003	AB-1, 50m x 0.2mm x 0.33um	Enquire
0132-5001	AB-1, 50m x 0.32mm x 0.17um	Enquire
0132-5005	AB-1, 50m x 0.32mm x 0.52um	Enquire
0132-5010	AB-1, 50m x 0.32mm x 1.05um	Enquire
0153-5030	AB-1, 50m x 0.53mm x 3.00um	Enquire
0153-5050	AB-1, 50m x 0.53mm x 5.00um	Enquire
0153-0501	AB-1, 5m x 0.53mm x 0.15um	Enquire
0153-0508	AB-1, 5m x 0.53mm x 0.88um	Enquire

0153-0526	AB-1, 5m x 0.53mm x 2.65um	Enquire
0153-0550	AB-1, 5m x 0.53mm x 5.0um	Enquire
0125-6002	AB-1, 60m x 0.25mm x 0.25um	Enquire
0125-6005	AB-1, 60m x 0.25mm x 0.50um	Enquire
0125-6010	AB-1, 60m x 0.25mm x 1.00um	Enquire
0132-6001	AB-1, 60m x 0.32mm x 0.10um	Enquire
0132-6002	AB-1, 60m x 0.32mm x 0.25um	Enquire
0132-6005	AB-1, 60m x 0.32mm x 0.50um	Enquire
0132-6010	AB-1, 60m x 0.32mm x 1.00um	Enquire
0132-6030	AB-1, 60m x 0.32mm x 3.00um	Enquire
0132-6050	AB-1, 60m x 0.32mm x 5.00um	Enquire
0153-6010	AB-1, 60m x 0.53mm x 1.00um	Enquire
0153-6030	AB-1, 60m x 0.53mm x 3.00um	Enquire
0153-6050	AB-1, 60m x 0.53mm x 5.00um	Enquire
0153-0715	AB-1, 7.5m x 0.53mm x 1.50um	Enquire
0153-0750	AB-1, 7.5m x 0.53mm x 5.00um	Enquire

### AB-1MS (100% Dimethylpolysiloxane)

#### Non Polar

Equivalent

to: BP-1, CP-Sil 5CB, DB-1, DB-1MS HP-1, HP-1MS  
Ultra-1, Rtx-1, SPB-1, ZB-1, ZB-1MS

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
1120-1503	AB-1MS, 15m x 0.20mm x 0.33um	Enquire
1125-1502	AB-1MS, 15m x 0.25mm x 0.25um	Enquire
1132-1502	AB-1MS, 15m x 0.32mm x 0.25um	Enquire
1120-2503	AB-1MS, 25m x 0.20mm x 0.33um	Enquire
1132-2505	AB-1MS, 25m x 0.32mm x 0.52um	Enquire
1125-3001	AB-1MS, 30m x 0.25mm x 0.10um	Enquire
1125-3002	AB-1MS, 30m x 0.25mm x 0.25um	Enquire
1125-3005	AB-1MS, 30m x 0.25mm x 0.50um	Enquire
1125-3010	AB-1MS, 30m x 0.25mm x 1.00um	Enquire
1132-3001	AB-1MS, 30m x 0.32mm x 0.10um	Enquire
1132-3002	AB-1MS, 30m x 0.32mm x 0.25um	Enquire
1132-3010	AB-1MS, 30m x 0.32mm x 1.00um	Enquire
1125-6002	AB-1MS, 60m x 0.25mm x 0.25um	Enquire
1132-6002	AB-1MS, 60m x 0.32mm x 0.25um	Enquire

### AB-1MS Builtin-Guard (100% Dimethylpolysiloxane)

#### Non Polar

Equivalent

to: BP-1, CP-Sil 5CB, DB-1, DB-1MS HP-1, HP-1MS  
Ultra-1, Rtx-1, SPB-1, ZB-1, ZB-1MS

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
1125-3002-G5	AB-1MS 30m,0.25mm,0.25um (5m Guard)	Enquire
1125-3002-G10	AB-1MS 30m,0.25mm,0.25um (10m Guard)	Enquire
1125-3005-G5	AB-1MS 30m,0.25mm,0.50um (5m Guard)	Enquire
1125-6002-G5	AB-1MS 60m,0.25mm,0.25um (5m Guard)	Enquire
1132-3002-G5	AB-1MS 30m,0.32mm,0.25um (5m Guard)	Enquire
1132-3005-G5	AB-1MS 30m,0.32mm,0.50um (5m Guard)	Enquire

## Ordering Information

### AB-PONA (100% Dimethylpolysiloxane)

#### Non Polar

Equivalent

to: DB-Petro, Petrocol DH, HP-PONA, SPB-1, 007-1,  
RTX-1, MXT-1

#### Part

#### Number

#### Description

#### UK

#### List

9006-PONA	AB-PONA 100m x 0.25mm x 0.5um	Enquire
9002-PONA	AB-PONA 50m x 0.20mm x 0.5um	Enquire

0532-5001	AB-5, 50m x 0.32mm x 0.17um	Enquire
0532-5005	AB-5, 50m x 0.32mm x 0.52um	Enquire
0532-5010	AB-5, 50m x 0.32mm x 1.05um	Enquire
0525-6001	AB-5, 60m x 0.25mm x 0.11um	Enquire
0525-6002	AB-5, 60m x 0.25mm x 0.25um	Enquire
0525-6005	AB-5, 60m x 0.25mm x 0.50um	Enquire
0525-6010	AB-5, 60m x 0.25mm x 1.00um	Enquire
0532-6001	AB-5, 60m x 0.32mm x 0.17um	Enquire
0532-6002	AB-5, 60m x 0.32mm x 0.25um	Enquire
0532-6010	AB-5, 60m x 0.32mm x 1.00um	Enquire
0553-6030	AB-5, 60m x 0.53mm x 3.00um	Enquire
0553-6050	AB-5, 60m x 0.53mm x 5.00um	Enquire

### AB-5 (5% Diphenyl 95% Dimethylpolysiloxane)

#### Non Polar

Equivalent

to: BP-5, CP-Sil 8CB, DB-5, HP-5, Ultra-2, Rtx-5,  
SPB-5, ZB-5

#### Part

#### Number

#### Description

#### UK

#### List

0553-1026	AB-5, 10m x 0.53mm x 2.65um	Enquire
0553-1050	AB-5, 10m x 0.53mm x 5.0um	Enquire
0520-1203	AB-5, 12m x 0.20mm x 0.33um	Enquire
0525-1501	AB-5, 15m x 0.25mm x 0.11um	Enquire
0525-1502	AB-5, 15m x 0.25mm x 0.25um	Enquire
0525-1505	AB-5, 15m x 0.25mm x 0.5um	Enquire
0525-1510	AB-5, 15m x 0.25mm x 1.0um	Enquire
0532-1501	AB-5, 15m x 0.32mm x 0.10um	Enquire
0532-1502	AB-5, 15m x 0.32mm x 0.25um	Enquire
0532-1505	AB-5, 15m x 0.32mm x 0.50um	Enquire
0532-1510	AB-5, 15m x 0.32mm x 1.00um	Enquire
0553-1515	AB-5, 15m x 0.53mm x 1.50um	Enquire
0553-1530	AB-5, 15m x 0.53mm x 3.00um	Enquire
0553-1550	AB-5, 15m x 0.53mm x 5.00um	Enquire
0520-1703	AB-5, 17m x 0.20mm x 0.33um	Enquire
0520-2501	AB-5, 25m x 0.20mm x 0.11um	Enquire
0520-2503	AB-5, 25m x 0.20mm x 0.33um	Enquire
0520-2505	AB-5, 25m x 0.20mm x 0.50um	Enquire
0532-2501	AB-5, 25m x 0.32mm x 0.17um	Enquire
0532-2505	AB-5, 25m x 0.32mm x 0.52um	Enquire
0532-2510	AB-5, 25m x 0.32mm x 1.05um	Enquire
0525-3001	AB-5, 30m x 0.25mm x 0.11um	Enquire
0525-3002	AB-5, 30m x 0.25mm x 0.25um	Enquire
0525-3005	AB-5, 30m x 0.25mm x 0.5um	Enquire
0525-3010	AB-5, 30m x 0.25mm x 1.00um	Enquire
0532-3001	AB-5, 30m x 0.32mm x 0.10um	Enquire
0532-3002	AB-5, 30m x 0.32mm x 0.25um	Enquire
0532-3005	AB-5, 30m x 0.32mm x 0.50um	Enquire
0532-3010	AB-5, 30m x 0.32mm x 1.00um	Enquire
0553-3002	AB-5, 30m x 0.53mm x 0.25um	Enquire
0553-3008	AB-5, 30m x 0.53mm x 0.88um	Enquire
0553-3015	AB-5, 30m x 0.53mm x 1.50um	Enquire
0553-3026	AB-5, 30m x 0.53mm x 2.65um	Enquire
0553-3030	AB-5, 30m x 0.53mm x 3.00um	Enquire
0553-3050	AB-5, 30m x 0.53mm x 5.00um	Enquire
0520-5001	AB-5, 50m x 0.20mm x 0.11um	Enquire
0520-5003	AB-5, 50m x 0.20mm x 0.33um	Enquire
0520-5005	AB-5, 50m x 0.20mm x 0.50um	Enquire

### AB-5MS (5% Diphenyl 95% Dimethylpolysiloxane)

#### Non Polar

Equivalent

to: DB-5MS, HP-5MS, Ultra-2, Rtx-5Sil MS, Rxi-5Sil MS,  
VF-5ms, CP-Sil 8 CB, BPX-5, ZB-5ms

#### Part

#### Number

#### Description

#### UK

#### List

1520-1203	AB-5MS, 12m x 0.20mm x 0.33um	Enquire
1525-1501	AB-5MS, 15m x 0.25mm x 0.10um	Enquire
1525-1502	AB-5MS, 15m x 0.25mm x 0.25um	Enquire
1525-1505	AB-5MS, 15m x 0.25mm x 0.50um	Enquire
1525-1510	AB-5MS, 15m x 0.25mm x 1.00um	Enquire
1532-1501	AB-5MS, 15m x 0.32mm x 0.10um	Enquire
1532-1502	AB-5MS, 15m x 0.32mm x 0.25um	Enquire
1532-1505	AB-5MS, 15m x 0.32mm x 0.50um	Enquire
1532-1510	AB-5MS, 15m x 0.32mm x 1.00um	Enquire
1520-1703	AB-5MS, 17m x 0.20mm x 0.33um	Enquire
1520-2503	AB-5MS, 25m x 0.20mm x 0.33um	Enquire
1532-2505	AB-5MS, 25m x 0.32mm x 0.52um	Enquire
1525-3001	AB-5MS, 30m x 0.25mm x 0.10um	Enquire
1525-3002	AB-5MS, 30m x 0.25mm x 0.25um	Enquire
1525-3005	AB-5MS, 30m x 0.25mm x 0.50um	Enquire
1525-3010	AB-5MS, 30m x 0.25mm x 1.00um	Enquire
1532-3001	AB-5MS, 30m x 0.32mm x 0.10um	Enquire
1532-3002	AB-5MS, 30m x 0.32mm x 0.25um	Enquire
1532-3005	AB-5MS, 30m x 0.32mm x 0.50um	Enquire
1532-3010	AB-5MS, 30m x 0.32mm x 1.00um	Enquire
1553-3005	AB-5MS, 30m x 0.53mm x 0.50um	Enquire
1553-3008	AB-5MS, 30m x 0.53mm x 0.88um	Enquire
1553-3015	AB-5MS, 30m x 0.53mm x 1.50um	Enquire
1520-5003	AB-5MS, 50m x 0.20mm x 0.33um	Enquire
1525-6001	AB-5MS, 60m x 0.25mm x 0.10um	Enquire
1525-6002	AB-5MS, 60m x 0.25mm x 0.25um	Enquire
1525-6005	AB-5MS, 60m x 0.25mm x 0.50um	Enquire
1525-6010	AB-5MS, 60m x 0.25mm x 1.00um	Enquire
1532-6001	AB-5MS, 60m x 0.32mm x 0.10um	Enquire
1532-6002	AB-5MS, 60m x 0.32mm x 0.25um	Enquire
1532-6005	AB-5MS, 60m x 0.32mm x 0.50um	Enquire
1532-6010	AB-5MS, 60m x 0.32mm x 1.00um	Enquire

## Ordering Information

### **AB-5MS Builtin-Guard (5% Diphenyl 95% Dimethylpolysiloxane)**

#### **Non Polar**

Equivalent

to: DB-5MS, HP-5MS, Ultra-2, Rtx-5Sil MS, Rxi-5Sil MS,  
VF-5ms, CP-Sil 8 CB, BPX-5, ZB-5ms

<b><u>Part Number</u></b>	<b><u>Description</u></b>	<b><u>UK List</u></b>
1525-3002-G5	AB-5MS 30m,0.25mm,0.25um (5m Guard)	Enquire
1525-3002-G10	AB-5MS 30m,0.25mm,0.25um (10m Guard)	Enquire
1525-3005-G5	AB-5MS 30m,0.25mm,0.50um (5m Guard)	Enquire
1525-6002-G5	AB-5MS 60m,0.25mm,0.25um (5m Guard)	Enquire
1532-3002-G5	AB-5MS 30m,0.32mm,0.25um (5m Guard)	Enquire
1532-3005-G5	AB-5MS 30m,0.32mm,0.50um (5m Guard)	Enquire

### **AB-35 (35% Phenyl 65% Dimethylpolysiloxane)**

#### **Mid Polar**

Equivalent

to: HP-35MS, DB-35MS, ZB-35, BPX-35, AT-35, Rtx-35,  
SPB-35

<b><u>Part Number</u></b>	<b><u>Description</u></b>	<b><u>UK List</u></b>
3525-15001	AB-35, 15m x 0.25mm x 0.10um	Enquire
3525-1502	AB-35, 15m x 0.25mm x 0.25um	Enquire
3525-1505	AB-35, 15m x 0.25mm x 0.50um	Enquire
3525-1510	AB-35, 15m x 0.25mm x 1.00um	Enquire
3532-1501	AB-35, 15m x 0.32mm x 0.10um	Enquire
3532-1502	AB-35, 15m x 0.32mm x 0.25um	Enquire
3532-1505	AB-35, 15m x 0.32mm x 0.50um	Enquire
3532-1510	AB-35, 15m x 0.32mm x 1.00um	Enquire
3532-1515	AB-35, 15m x 0.32mm x 1.50um	Enquire
3553-1501	AB-35, 15m x 0.53mm x 0.10um	Enquire
3553-1502	AB-35, 15m x 0.53mm x 0.25um	Enquire
3553-1505	AB-35, 15m x 0.53mm x 0.50um	Enquire
3553-1510	AB-35, 15m x 0.53mm x 1.00um	Enquire
3553-1515	AB-35, 15m x 0.53mm x 1.50um	Enquire
3525-3001	AB-35, 30m x 0.25mm x 0.15um	Enquire
3525-3002	AB-35, 30m x 0.25mm x 0.25um	Enquire
3525-3005	AB-35, 30m x 0.25mm x 0.50um	Enquire
3525-3010	AB-35, 30m x 0.25mm x 1.00um	Enquire
3532-3001	AB-35, 30m x 0.32mm x 0.10um	Enquire
3532-3002	AB-35, 30m x 0.32mm x 0.25um	Enquire
3532-3005	AB-35, 30m x 0.32mm x 0.50um	Enquire
3532-3010	AB-35, 30m x 0.32mm x 1.00um	Enquire
3532-3015	AB-35, 30m x 0.32mm x 1.50um	Enquire
3553-3001	AB-35, 30m x 0.53mm x 0.10um	Enquire
3553-3002	AB-35, 30m x 0.53mm x 0.25um	Enquire
3553-3005	AB-35, 30m x 0.53mm x 0.50um	Enquire
3553-3010	AB-35, 30m x 0.53mm x 1.00um	Enquire
3553-3015	AB-35, 30m x 0.53mm x 1.50um	Enquire
3525-6001	AB-35, 60m x 0.25mm x 0.10um	Enquire
3525-6002	AB-35, 60m x 0.25mm x 0.25um	Enquire
3525-6005	AB-35, 60m x 0.25mm x 0.50um	Enquire
3525-6010	AB-35, 60m x 0.25mm x 1.00um	Enquire
3532-6001	AB-35, 60m x 0.32mm x 0.10um	Enquire
3532-6002	AB-35, 60m x 0.32mm x 0.25um	Enquire
3532-6005	AB-35, 60m x 0.32mm x 0.50um	Enquire

3532-6010	AB-35, 60m x 0.32mm x 1.00um	Enquire
3532-6015	AB-35, 60m x 0.32mm x 1.50um	Enquire
3553-6001	AB-35, 60m x 0.53mm x 0.10um	Enquire
3553-6002	AB-35, 60m x 0.53mm x 0.25um	Enquire
3553-6005	AB-35, 60m x 0.53mm x 0.50um	Enquire
3553-6010	AB-35, 60m x 0.53mm x 1.00um	Enquire
3553-6015	AB-35, 60m x 0.53mm x 1.50um	Enquire

### **AB-624 (6% Cyanopropylphenyl 94% Dimethylpolysiloxane)**

#### **Mid Polar**

Equivalent

to: DB-624, ZB-624, AT-624, Rtx-624

<b><u>Part Number</u></b>	<b><u>Description</u></b>	<b><u>UK List</u></b>
6225-3014	AB-624, 30m x 0.25mm x 1.40um	Enquire
6232-3018	AB-624, 30m x 0.32mm x 1.80um	Enquire
6253-3030	AB-624, 30m x 0.53mm x 3.00um	Enquire
6225-6014	AB-624, 60m x 0.25mm x 1.40um	Enquire
6232-6018	AB-624, 60m x 0.32mm x 1.80um	Enquire

### **AB-1301 (6% Cyanopropylphenyl 94% Dimethylpolysiloxane)**

#### **Mid Polar**

Equivalent

to: HP-1301, DB-1301, Rtx-1301

<b><u>Part Number</u></b>	<b><u>Description</u></b>	<b><u>UK List</u></b>
6025-1501	AB-1301, 15m x 0.25mm x 0.10um	Enquire
6025-1502	AB-1301, 15m x 0.25mm x 0.25um	Enquire
6025-1505	AB-1301, 15m x 0.25mm x 0.50um	Enquire
6025-1510	AB-1301, 15m x 0.25mm x 1.00um	Enquire
6032-1501	AB-1301, 15m x 0.32mm x 0.10um	Enquire
6032-1502	AB-1301, 15m x 0.32mm x 0.25um	Enquire
6032-1505	AB-1301, 15m x 0.32mm x 0.50um	Enquire
6032-1510	AB-1301, 15m x 0.32mm x 1.00um	Enquire
6053-1501	AB-1301, 15m x 0.53mm x 0.10um	Enquire
6053-1502	AB-1301, 15m x 0.53mm x 0.25um	Enquire
6053-1505	AB-1301, 15m x 0.53mm x 0.50um	Enquire
6053-1510	AB-1301, 15m x 0.53mm x 1.00um	Enquire
6053-1515	AB-1301, 15m x 0.53mm x 1.50um	Enquire
6053-1530	AB-1301, 15m x 0.53mm x 3.00um	Enquire
6025-3001	AB-1301, 30m x 0.25mm x 0.10um	Enquire
6025-3002	AB-1301, 30m x 0.25mm x 0.25um	Enquire
6025-3005	AB-1301, 30m x 0.25mm x 0.50um	Enquire
6025-3010	AB-1301, 30m x 0.25mm x 1.00um	Enquire
6032-3001	AB-1301, 30m x 0.32mm x 0.10um	Enquire
6032-3002	AB-1301, 30m x 0.32mm x 0.25um	Enquire
6032-3005	AB-1301, 30m x 0.32mm x 0.50um	Enquire
6032-3010	AB-1301, 30m x 0.32mm x 1.00um	Enquire
6032-3018	AB-1301, 30m x 0.32mm x 1.80um	Enquire
6053-3001	AB-1301, 30m x 0.53mm x 0.10um	Enquire
6053-3002	AB-1301, 30m x 0.53mm x 0.25um	Enquire
6053-3005	AB-1301, 30m x 0.53mm x 0.50um	Enquire
6053-3010	AB-1301, 30m x 0.53mm x 1.00um	Enquire

## Ordering Information

6053-3015	AB-1301, 30m x 0.53mm x 1.50um	Enquire	6125-6005	AB-1701, 60m x 0.25mm x 0.5um	Enquire
6053-3030	AB-1301, 30m x 0.53mm x 3.00um	Enquire	6125-6010	AB-1701, 60m x 0.25mm x 1.00um	Enquire
6025-6001	AB-1301, 60m x 0.25mm x 0.10um	Enquire	6132-6001	AB-1701, 60m x 0.32mm x 0.10um	Enquire
6025-6002	AB-1301, 60m x 0.25mm x 0.25um	Enquire	6132-60015	AB-1701, 60m x 0.32mm x 0.15um	Enquire
6025-6005	AB-1301, 60m x 0.25mm x 0.50um	Enquire	6132-6002	AB-1701, 60m x 0.32mm x 0.25um	Enquire
6025-6010	AB-1301, 60m x 0.25mm x 1.00um	Enquire	6132-6005	AB-1701, 60m x 0.32mm x 0.50um	Enquire
6032-6001	AB-1301, 60m x 0.32mm x 0.10um	Enquire	6132-6010	AB-1701, 60m x 0.32mm x 1.00um	Enquire
6032-6002	AB-1301, 60m x 0.32mm x 0.25um	Enquire	6153-6001	AB-1701, 60m x 0.53mm x 0.1um	Enquire
6032-6005	AB-1301, 60m x 0.32mm x 0.50um	Enquire	6153-6002	AB-1701, 60m x 0.53mm x 0.25um	Enquire
6032-6010	AB-1301, 60m x 0.32mm x 1.00um	Enquire	6153-6005	AB-1701, 60m x 0.53mm x 0.5um	Enquire
6053-6001	AB-1301, 60m x 0.53mm x 0.10um	Enquire	6153-6010	AB-1701, 60m x 0.53mm x 1.00um	Enquire
6053-6002	AB-1301, 60m x 0.53mm x 0.25um	Enquire	6153-6015	AB-1701, 60m x 0.53mm x 1.50um	Enquire
6053-6005	AB-1301, 60m x 0.53mm x 0.50um	Enquire			
6053-6010	AB-1301, 60m x 0.53mm x 1.00um	Enquire			
6053-6015	AB-1301, 60m x 0.53mm x 1.50um	Enquire			
6053-6030	AB-1301, 60m x 0.53mm x 3.00um	Enquire			

### AB-1701 (14% Cyanopropylphenyl 86% Dimethylpolysiloxane)

#### Mid Polar

Equivalent to: HP-1701, DB-1701, ZB-1701, BP-10, CP-Sil 19CB

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
6125-1501	AB-1701, 15m x 0.25mm x 0.10um	Enquire
6125-15015	AB-1701, 15m x 0.25mm x 0.15um	Enquire
6125-1502	AB-1701, 15m x 0.25mm x 0.25um	Enquire
6125-1505	AB-1701, 15m x 0.25mm x 0.50um	Enquire
6125-1510	AB-1701, 15m x 0.25mm x 1.00um	Enquire
6132-1501	AB-1701, 15m x 0.32mm x 0.10um	Enquire
6132-15015	AB-1701, 15m x 0.32mm x 0.15um	Enquire
6132-1502	AB-1701, 15m x 0.32mm x 0.25um	Enquire
6132-1505	AB-1701, 15m x 0.32mm x 0.50um	Enquire
6132-1510	AB-1701, 15m x 0.32mm x 1.00um	Enquire
6153-1501	AB-1701, 15m x 0.53mm x 0.10um	Enquire
6153-1502	AB-1701, 15m x 0.53mm x 0.25um	Enquire
6153-1505	AB-1701, 15m x 0.53mm x 0.50um	Enquire
6153-1510	AB-1701, 15m x 0.53mm x 1.00um	Enquire
6153-1515	AB-1701, 15m x 0.53mm x 1.50um	Enquire
6125-3001	AB-1701, 30m x 0.25mm x 0.10um	Enquire
6125-30015	AB-1701, 30m x 0.25mm x 0.15um	Enquire
6125-3002	AB-1701, 30m x 0.25mm x 0.25um	Enquire
6125-3005	AB-1701, 30m x 0.25mm x 0.50um	Enquire
6125-3010	AB-1701, 30m x 0.25mm x 1.00um	Enquire
6132-3001	AB-1701, 30m x 0.32mm x 0.10um	Enquire
6132-30015	AB-1701, 30m x 0.32mm x 0.15um	Enquire
6132-3002	AB-1701, 30m x 0.32mm x 0.25um	Enquire
6132-3005	AB-1701, 30m x 0.32mm x 0.50um	Enquire
6132-3010	AB-1701, 30m x 0.32mm x 1.00um	Enquire
6132-3018	AB-1701, 30m x 0.32mm x 1.80um	Enquire
6153-3001	AB-1701, 30m x 0.53mm x 0.10um	Enquire
6153-3002	AB-1701, 30m x 0.53mm x 0.25um	Enquire
6153-3005	AB-1701, 30m x 0.53mm x 0.50um	Enquire
6153-3010	AB-1701, 30m x 0.53mm x 1.00um	Enquire
6153-3015	AB-1701, 30m x 0.53mm x 1.50um	Enquire
6125-6001	AB-1701, 60m x 0.25mm x 0.10um	Enquire
6125-60015	AB-1701, 60m x 0.25mm x 0.15um	Enquire
6125-6002	AB-1701, 60m x 0.25mm x 0.25um	Enquire

### AB-225 (50% Cyanopropylphenyl - Dimethylpolysiloxane)

#### Mid-High Polarity

Equivalent to: SP-2330, CP Sil 43 CB, Rtx-225, BP-225, OV-225, 007-225, AT-225

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
2525-1501	AB-225, 15m x 0.25mm x 0.10	Enquire
2525-15015	AB-225, 15m x 0.25mm x 0.15	Enquire
2525-1502	AB-225, 15m x 0.25mm x 0.25	Enquire
2525-1505	AB-225, 15m x 0.25mm x 0.50	Enquire
2532-1501	AB-225, 15m x 0.32mm x 0.10	Enquire
2532-15015	AB-225, 15m x 0.32mm x 0.15	Enquire
2532-1502	AB-225, 15m x 0.32mm x 0.25	Enquire
2532-1505	AB-225, 15m x 0.32mm x 0.50	Enquire
2553-1501	AB-225, 15m x 0.53mm x 0.10	Enquire
2553-1502	AB-225, 15m x 0.53mm x 0.25	Enquire
2553-1505	AB-225, 15m x 0.53mm x 0.50	Enquire
2553-1510	AB-225, 15m x 0.53mm x 1.00	Enquire
2525-3001	AB-225, 30m x 0.25mm x 0.10	Enquire
2525-30015	AB-225, 30m x 0.25mm x 0.15	Enquire
2525-3002	AB-225, 30m x 0.25mm x 0.25	Enquire
2525-3005	AB-225, 30m x 0.25mm x 0.50	Enquire
2532-30015	AB-225, 30m x 0.32mm x 0.15	Enquire
2532-3002	AB-225, 30m x 0.32mm x 0.25	Enquire
2532-3005	AB-225, 30m x 0.32mm x 0.50	Enquire
2553-3001	AB-225, 30m x 0.53mm x 0.10	Enquire
2553-3002	AB-225, 30m x 0.53mm x 0.25	Enquire
2553-3005	AB-225, 30m x 0.53mm x 0.50	Enquire
2553-3010	AB-225, 30m x 0.53mm x 1.00	Enquire

## Ordering Information

### AB-CarboWax 20M (Polyethylene Glycol)

**Polar**

Equivalent to: HP-20M, DB-CAM

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
2253-1013	AB-CarboWax 10m x 0.53mm x 1.33um	Enquire
2220-2501	AB-CarboWax 25m x 0.20mm x 0.10um	Enquire
2232-2503	AB-CarboWax 25m x 0.32mm x 0.30um	Enquire
2225-3002	AB-CarboWax 30m x 0.25mm x 0.25um	Enquire
2253-3013	AB-CarboWax 30m x 0.53mm x 1.33um	Enquire
2220-5001	AB-CarboWax 50m x 0.20mm x 0.10um	Enquire
2232-5003	AB-CarboWax 50m x 0.32mm x 0.30um	Enquire

### AB-Inowax (Polyethylene Glycol)

**Polar**

Equivalent to: HP-Innowax, DB-WAX, DB-WAXetr Supelco WAX,

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
2025-1502	AB-INOWAX, 15m x 0.25mm x 0.25um	Enquire
2025-1505	AB-INOWAX, 15m x 0.25mm x 0.50um	Enquire
2032-1502	AB-INOWAX, 15m x 0.32mm x 0.25um	Enquire
2032-1505	AB-INOWAX, 15m x 0.32mm x 0.50um	Enquire
2053-1510	AB-INOWAX, 15m x 0.53mm x 1.00um	Enquire
2020-2502	AB-INOWAX, 25m x 0.20mm x 0.20um	Enquire
2020-2504	AB-INOWAX, 25m x 0.20mm x 0.40um	Enquire
2025-3001	AB-INOWAX, 30m x 0.25mm x 0.15um	Enquire
2025-3002	AB-INOWAX, 30m x 0.25mm x 0.25um	Enquire
2025-3005	AB-INOWAX, 30m x 0.25mm x 0.50um	Enquire
2032-3001	AB-INOWAX, 30m x 0.32mm x 0.15um	Enquire
2032-3002	AB-INOWAX, 30m x 0.32mm x 0.25um	Enquire
2032-3005	AB-INOWAX, 30m x 0.32mm x 0.50um	Enquire
2053-3010	AB-INOWAX, 30m x 0.53mm x 1.00um	Enquire
2020-5002	AB-INOWAX, 50m x 0.20mm x 0.20um	Enquire
2020-5004	AB-INOWAX, 50m x 0.20mm x 0.40um	Enquire
2025-6001	AB-INOWAX, 60m x 0.25mm x 0.15um	Enquire
2025-6002	AB-INOWAX, 60m x 0.25mm x 0.25um	Enquire
2025-6005	AB-INOWAX, 60m x 0.25mm x 0.50um	Enquire
2032-6001	AB-INOWAX, 60m x 0.32mm x 0.15um	Enquire
2032-6002	AB-INOWAX, 60m x 0.32mm x 0.25um	Enquire
2032-6005	AB-INOWAX, 60m x 0.32mm x 0.50um	Enquire
2053-6010	AB-INOWAX, 60m x 0.53mm x 1.00um	Enquire

### AB-FFAP (Polyethylene Glycol Acid Modified)

**Polar**

Equivalent to: HP-FFAP, DB-FFAP, CP-WAX 58 CB, BP21, AT-100

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
2153-1010	AB-FFAP, 10m x 0.53mm x 1.00um	Enquire
2125-1502	AB-FFAP, 15m x 0.25mm x 0.25um	Enquire
2132-1502	AB-FFAP, 15m x 0.32mm x 0.25um	Enquire
2153-1510	AB-FFAP, 15m x 0.53mm x 1.00um	Enquire
2120-2503	AB-FFAP, 25m x 0.20mm x 0.30um	Enquire
2132-2505	AB-FFAP, 25m x 0.32mm x 0.50um	Enquire
2125-3002	AB-FFAP, 30m x 0.25mm x 0.25um	Enquire
2132-3002	AB-FFAP, 30m x 0.32mm x 0.25um	Enquire
2132-3005	AB-FFAP, 30m x 0.32mm x 0.50um	Enquire
2153-3010	AB-FFAP, 30m x 0.53mm x 1.00um	Enquire
2120-5003	AB-FFAP, 50m x 0.20mm x 0.30um	Enquire
2132-5005	AB-FFAP, 50m x 0.32mm x 0.50um	Enquire

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
8132-1508	AB-PLOT Al2O3 KCl 15m x 0.32mm x 8.00um	Enquire
8153-1515	AB-PLOT Al2O3 KCl 15m x 0.53mm x 15.0um	Enquire
8125-2505	AB-PLOT Al2O3 KCl 25m x 0.25mm x 5.00um	Enquire
8125-3005	AB-PLOT Al2O3 KCl 30m x 0.25mm x 5.00um	Enquire
8132-3008	AB-PLOT Al2O3 KCl 30m x 0.32mm x 8.00um	Enquire
8153-3015	AB-PLOT Al2O3 KCl 30m x 0.53mm x 15.0um	Enquire
8132-5008	AB-PLOT Al2O3 KCl 50m x 0.32mm x 8.00um	Enquire
8153-5015	AB-PLOT Al2O3 KCl 50m x 0.53mm x 15.0um	Enquire
8153-6015	AB-PLOT Al2O3 KCl 60m x 0.53mm x 15.0um	Enquire

### AB-PLOT Al2O3-KCl

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
8332-1508	AB-PLOT Al2O3 "M", 15m x 0.32mm x 8.00um	Enquire
8353-1515	AB-PLOT Al2O3 "M", 15m x 0.53mm x 15.0um	Enquire
8325-2505	AB-PLOT Al2O3 "M", 25m x 0.25mm x 5.00um	Enquire
8325-3005	AB-PLOT Al2O3 "M", 30m x 0.25mm x 5.00um	Enquire
8332-3008	AB-PLOT Al2O3 "M", 30m x 0.32mm x 8.00um	Enquire
8353-3015	AB-PLOT Al2O3 "M", 30m x 0.53mm x 15.0um	Enquire
8332-5008	AB-PLOT Al2O3 "M", 50m x 0.32mm x 8.00um	Enquire
8353-5015	AB-PLOT Al2O3 "M", 50m x 0.53mm x 15.0um	Enquire
8353-6015	AB-PLOT Al2O3 "M", 60m x 0.53mm x 15.0um	Enquire

### AB-PLOT Al2O3-M

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
8232-1508	AB-PLOT Al2O3 "S", 15m x 0.32mm x 8.00um	Enquire
8253-1515	AB-PLOT Al2O3 "S", 15m x 0.53mm x 15.0um	Enquire
8225-2505	AB-PLOT Al2O3 "S", 25m x 0.25mm x 5.00um	Enquire
8225-3005	AB-PLOT Al2O3 "S", 30m x 0.25mm x 5.00um	Enquire
8232-3008	AB-PLOT Al2O3 "S", 30m x 0.32mm x 8.00um	Enquire
8253-3015	AB-PLOT Al2O3 "S", 30m x 0.53mm x 15.0um	Enquire
8232-5008	AB-PLOT Al2O3 "S", 50m x 0.32mm x 8.00um	Enquire

### AB-PLOT Al2O3-S

<u>Part Number</u>	<u>Description</u>	<u>UK List</u>
8232-1508	AB-PLOT Al2O3 "S", 15m x 0.32mm x 8.00um	Enquire
8253-1515	AB-PLOT Al2O3 "S", 15m x 0.53mm x 15.0um	Enquire
8225-2505	AB-PLOT Al2O3 "S", 25m x 0.25mm x 5.00um	Enquire
8225-3005	AB-PLOT Al2O3 "S", 30m x 0.25mm x 5.00um	Enquire
8232-3008	AB-PLOT Al2O3 "S", 30m x 0.32mm x 8.00um	Enquire
8253-3015	AB-PLOT Al2O3 "S", 30m x 0.53mm x 15.0um	Enquire
8232-5008	AB-PLOT Al2O3 "S", 50m x 0.32mm x 8.00um	Enquire

## Ordering Information

8253-5015	AB-PLOT Al <sub>2</sub> O <sub>3</sub> "S", 50m x 0.53mm x 15.0um	Enquire
8253-6015	AB-PLOT Al <sub>2</sub> O <sub>3</sub> "S", 60m x 0.53mm x 15.0um	Enquire

### **AB-PLOT Molsieve**

Equivalent to: HP-PLOT Molsieve, CP-PLOT Molsieve

<b><u>Part Number</u></b>	<b><u>Description</u></b>	<b><u>UK List</u></b>
8432-1512	AB-PLOT MoleSieve, 15m x 0.32mm x 12.0um	Enquire
8432-1525	AB-PLOT MoleSieve, 15m x 0.32mm x 25.0um	Enquire
8453-1525	AB-PLOT MoleSieve, 15m x 0.53mm x 25.0um	Enquire
8453-1550	AB-PLOT MoleSieve, 15m x 0.53mm x 50.0um	Enquire
8432-3012	AB-PLOT MoleSieve, 30m x 0.32mm x 12.0um	Enquire
8432-3025	AB-PLOT MoleSieve, 30m x 0.32mm x 25.0um	Enquire
8453-3025	AB-PLOT MoleSieve, 30m x 0.53mm x 25.0um	Enquire
8453-3050	AB-PLOT MoleSieve, 30m x 0.53mm x 50.0um	Enquire
8432-5025	AB-PLOT MoleSieve, 50m x 0.32mm x 25.0um	Enquire
8453-5050	AB-PLOT MoleSieve, 50m x 0.53mm x 50.0um	Enquire
8432-0512	AB-PLOT MoleSieve, 5m x 0.32mm x 12.0um	Enquire
8453-0525	AB-PLOT MoleSieve, 5m x 0.53mm x 25.0um	Enquire
8453-6025	AB-PLOT MoleSieve, 60m x 0.32mm x 25.0um	Enquire
8453-6050	AB-PLOT MoleSieve, 60m x 0.53mm x 50.0um	Enquire

### **AB-PLOT-Q**

Equivalent to: HP-PLOT-Q, GS-Q, CP-PoraPLOT Q, Rt-Q

<b><u>Part Number</u></b>	<b><u>Description</u></b>	<b><u>UK List</u></b>
8632-1515	AB-PLOT Q, 15m x 0.32mm x 15.0um	Enquire
8653-1530	AB-PLOT Q, 15m x 0.53mm x 30.0um	Enquire
8632-3015	AB-PLOT Q, 30m x 0.32mm x 15.0um	Enquire
8653-3030	AB-PLOT Q, 30m x 0.53mm x 30.0um	Enquire
8653-6030	AB-PLOT Q, 60m x 0.53mm x 30.0um	Enquire

### **AB-PLOT-U**

Equivalent to: HP-PLOT U, CP-PoraPLOT U

<b><u>Part Number</u></b>	<b><u>Description</u></b>	<b><u>UK List</u></b>
8732-1510	AB-PLOT U, 15m x 0.32mm x 10.0um	Enquire
8753-1520	AB-PLOT U, 15m x 0.53mm x 20.0um	Enquire
8732-3010	AB-PLOT U, 30m x 0.32mm x 10.0um	Enquire
8732-3020	AB-PLOT U, 30m x 0.32mm x 20.0um	Enquire
8753-3020	AB-PLOT U, 30m x 0.53mm x 20.0um	Enquire